Glucosefine®

KONTROLL-LÖSUNG

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CONTROL SOLUTION

Intended use

The CONTROL SOLUTION is used with the Glucosefine® Test Strips. This ensures that the meter and strips are working together properly and/or the test is being performed correctly for the measurement of glucose in your blood. It is very important that you test with CONTROL SOLUTION routinely to make sure your results are accurate. The control solutions are for testing outside the body (in vitro diagnostic use only). Glucosefine® system is for self-testing.

Quality Control

Acceptable performance is assured when the control test result falls within the appropriate range printed on the test strip vial.

Summary

The CONTROL SOLUTION contains a known amount of glucose that reacts with the test strip. A test with CONTROL SOLUTION is similar to that of a blood test. All you do is use your CONTROL SOLUTION instead of a drop of blood. The CONTROL SOLUTION result should fall within the appropriate range printed on the test strip vial.

NOTE: If the control solution test is not in range, do NOT continue with a blood glucose test. Please contact customer service of Metrado GmbH or the distributor in your country.

When to run a control test

- Before using a new box of strips.
- Whenever you feel your glucose levels are inaccurate.
- Before adjusting medication when a blood glucose test result is too low or too high.
- If you feel the meter or test strips have been compromised in any way.

Storage and handling

- Make sure the CONTROL SOLUTION bottle is tightly closed and store at room temperature between 8 - 30 °C.
 - If the CONTROL SOLUTION is hot or cold, do not use until the solution has remained at room temperature for at least 10 minutes. Perform test between 20 - 25 °C.
 - Use before the expiration date shown on the bottle.
- Once opened, use within 3 months of that time frame. Make sure to discard after three months.

Testing procedure

- Accurate results depend upon carefully following all instructions for use.
1. It is recommended that the meter, control solution, and test strips are at room temperature before testing. Test between 20 - 25 °C.
 2. Prior to use, shake the bottle vigorously.
 3. Discard first drop of CONTROL SOLUTION.
 4. Wipe the tip clean.
 5. Make sure the test strip is inserted fully into the test slot of the meter.
 6. The meter will automatically turn on and display the code number. Make sure the code number matches the code number on the test strip vial.
 7. Specify that it's a CONTROL SOLUTION test and not a regular blood glucose test by pressing the up or down button. A medication bottle will appear on the top of the screen.
 8. Place a drop of solution on a clean, dry surface (for example: the lid of your test strip vial). Pick up the meter and touch the end of the test strip to the drop. Be sure to hold the test strip to the drop until the "beep" sound.
 9. A result will appear in 5 seconds.
 10. Close the bottle immediately after use.

Reagents

Each Glucosefine® CONTROL SOLUTION is provided as a buffered, stabilized, aqueous solution of D-glucose containing preservative and other non-reactive ingredients.


Normal: 0.12 wt% glucose 99.88 wt% non-reactive ingredients
High: 0.20 wt% glucose 99.80 wt% non-reactive ingredients

Precautions

- The CONTROL SOLUTION is for In-Vitro Diagnostic use only. This means that it is only used for testing outside of the body.
- The CONTROL SOLUTION is made for the Glucosefine® system, do not use any other brand of control solution.
- Check the expiration date on the bottle. Do not use if expired.
- Once opened, only use for the next three months. Record the discard date on the control solution bottle. Discard after three months.
- Do not swallow. Not for human or animal consumption.

Expected Results

At room temperature, CONTROL SOLUTION test results fall within the expected range printed on the test strip vial. For example only: If you are using the "normal" CONTROL SOLUTION, according to this particular vial of strips, your meter should show a number between 88 and 132mg/dL (4,8 and 7,3 mmol/L).

If test results fall outside the expected range, repeat the test. Do not continue with your blood glucose test if the control solution test results are continuously out of range. Results that fall outside the expected range may indicate:

- Test not performed at room temperature.
- An error in the way you are doing the test.
- Expired or contaminated CONTROL SOLUTION.
- A problem with the meter.
- CONTROL SOLUTION that is too warm or cold.
- Failure to discard the first drop of CONTROL SOLUTION and to wipe the bottle tip clean.
- Failure to shake the CONTROL SOLUTION adequately.

NOTE: Please refer to the table below to identify symbols.

- Consult operating instructions
- Expiration date
- This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices.
- For in vitro diagnostic use
- Batch code
- Catalogue number
- Manufacturing date
- Self-testing Used
- Temperature limitation
- Attention, see instructions for use
- Manufacturer
- Discard 3 months after opening
- Disposed of with household waste. Please observe the local regulations as these may vary from country to country.

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METRADO GmbH