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Glucosefine® Glucose Test Strips

IMPORTANT: Please read this information and your Glucosefine® operation manual before using the Glucosefine® glucose test strips.

Intended Use:

The Glucosefine® Glucose Test Strips are used with Glucosefine® Glucose Test Meter for the measurement of glucose in your blood. The Glucosefine® Blood Glucose Monitoring System is plasma-calibrated to allow comparison of results with the laboratory method. The Glucosefine® Glucose Test Strips are for testing outside the body (in vitro diagnostic use only). They should be used only for testing blood glucose (sugar) and only with fresh capillary whole blood samples. Glucosefine® Blood Glucose Monitoring System is for self-testing by people with diabetes and by healthcare professionals. However, the self-testing should not take the place of overall diabetes management assessment from your doctor.

General:

Glucosefine® glucose test strips adopt a bio-sensor technology composed of electro-chemical components and require a small (0.5µl) volume of blood giving you less pain. Glucosefine® can store 365 readings in memory with date & time, you can easily store your readings on your PC (optional). This will provide an easier and better way to control your glucose levels.

Storage and Handling:

- Store the Glucosefine® Glucose Test Strip vials in a cool and dry place between 2-30°C (36-86°F). Keep out of direct sunlight. Do not freeze.
- Store Glucose Test Strips in their original vial only. Do not mix the Glucose Test Strips in new vials or in any other container.
- Immediately replace the cap and close tightly after removing a Glucosefine® Glucose Test Strip.
- Use each Glucose Test Strip immediately after removing it from the vial.
- Do not use Glucose Test Strips after the expiry date printed on the package or vial since it may cause inaccurate results.
- Make a notation of the discard date on the vial label when you first open it.
- Discard remaining Glucosefine® Glucose Test Strips 6 months after first opening the vial.
- Avoid getting dirt, food, and water on the Glucose Test Strip. Do not handle Glucose Test Strips with wet hands.
- Do not bend, cut, or alter any Glucosefine® Glucose Test Strip.
- Glucosefine® Glucose Test Strips are for single use only. Do not re-use.
- Do not perform tests at temperatures below 10°C (50°F) or above 40°C (104 °F). Do not perform tests with humidity below 10 % or above 90 %.
- Be careful not to pollute the color coding label. (Polluted color coding label due to dust and other substances)
- Discard the used Glucose Test Strip carefully, to prevent any infection.
- Keep the Test Strips away from children!

Warning! To avoid false readings, use only a Glucosefine® Meter to measure glucose with your Glucosefine® Glucose Test Strips.

Care Precautions:

- Do not reuse glucose test strips. Single use only.
- If you experience any symptoms that are not consistent with your blood glucose test results and you have followed the instructions described in your Glucosefine® Meter Manual, call your physician.
- Do not make significant changes to your diabetes control program without consulting your physician.
- Do not ignore physical symptoms without consulting your physician.

Test Principle:

Glucose in the blood sample will react to the electrodes in the glucose test strip, generating an electrical current that will stimulate a chemical reaction. This reaction is measured by the Glucosefine® meter and displayed as your blood glucose result.

Reagent Composition:

Each Glucosefine® glucose test strip contains:
Glucose Oxidase: 1.5 Units
Mediator: 100 µg

The procedure for Blood Glucose Measurement:

Materials needed: Glucosefine® Test Strips, Glucosefine® meter, operation manual, lancing device, sterile lancet

Obtaining Blood Sample:

The Glucosefine® Glucose Test Strips are designed to be used with fresh capillary whole blood. You may obtain a blood drop from either a fingertip or an alternative site (see meter manual).

To obtain a drop of blood, follow these steps :

Step 1: Cleanliness

Wash your hands and test site with warm water. Dry it thoroughly. (Dirt, oil, lotion, or sweat may affect the test result).

Step 2: Lancing

Prepare the lancing device. Insert a lancet into the lancing device. The device holds, positions, and controls how deeply the lancet goes into the skin. (Refer to the owner’s manual for more information).


Step 3: Insert the Glucose Test Strip

Remove a test strip from the vial. Recap the vial immediately to prevent moisture from affecting the other strips. Insert the test strip into the port of the meter, with arrow sign facing up and showing towards the meter. The meter will automatically turn on and display the code number. Make sure that the code on your meter display matches the code number printed on the test strip vial. If not matching repeat this step.

The test strip icon with a drop of blood will be flashing at the top of the meter display. You are now ready for testing.

Step 4: Obtain and apply the blood sample

Place the lancing device against the area to be pricked, for example the finger tip. For a correct test result, an amount of 0.5µl is necessary. While the symbol for apply is flashing, place the drop of blood on your finger against the slim end of the test strip until the window is completely filled. Do not press your finger against the test strip. Do not try to wipe off the blood with the test strip. A beep sound confirms, when the window is correctly filled. If you have a poor amount of blood: Run a new test whith a NEW test strip!

Step 5: Measurement result in 5 seconds

The meter automatically starts testing. It counts down from 5 to 1, displays the measurement result and stores it automatically (To store with User’s Activity and User-ID please refer to your meter manual). Use the test strip ejector to remove the test strip. The meter will automatically turn off.

Important Information About Using Alternate Site Testing:

- Alternate sites where you can test are dorsal hand, ventral palm, upper arm, forearm, calf, and thigh.
- Under certain conditions, blood glucose test results obtained using samples taken from your alternate sites may differ significantly from fingertip samples.
- The conditions in which these differences are more likely to occur are when your blood glucose is changing rapidly such as following a meal, as insulin dose, or associated with physical exercise.
- When blood glucose is changing rapidly, fingertip samples show these changes more quickly than alternate sites samples.
- When your blood glucose is falling, testing with a fingertip sample may identify a hypoglycemic (low blood sugar) level sooner than a test with a alternate sites sample.
- Use alternate sites samples only for testing prior to or more than two hours after meal, insulin doses, or physical exercise.
- Testing performed within two hours after meals, insulin doses, or physical exercise or whenever you feel that your glucose levels may be changing rapidly should be done from the fingertip.

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